

Glass DC Nebulizers

SeaSpray

400µL/min uptake A52-1-USS04

MicroMist

400µL/min uptake A52-1-UM04

200µL/min uptake A52-1-UM02

100µL/min uptake A52-1-UM01



DC Nebulizers for HF

OpalMist

50µL/min uptake A52-1-PFA005



DC Nebulizers for HE-SIS

MicroMist HE

200µL/min uptake C52-1-UFT02



Spray Chambers

High Efficiency Sample Introduction System

HE-SIS for Nu Vitesse KT-1219



Twister

50mL Cyclonic 20-809-0495HE



Torches

Quartz Torch

with 1.5mm Injector 30-807-0560



Viton O-rings for FDT/SDT ball joints (PKT 10)

70-V-011



Cones

Sampler Cones

Nickel Sampler Version A, Wet Plasma NU1011A-Ni

Nickel Dry Sampler with Copper Core NU1011B-Ni

Skimmer Cones

Nickel Skimmer Version A, Wet Plasma NU1004A-Ni

Nickel Dry Skimmer NU1014B-Ni



Common Accessories

Guardian In-Line Filters


Guardian DC in-line non-return gas filter 70-803-1942



Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



Magnifier Inspection Tool

70-803-1923

